

## Bp Premier Spectra SP1 Revision 1 Release Notes

---

These release notes describe the changes made to Bp Premier since Spectra SP1 build 1.13.1.1075.

Release Date	28 May 2025
Release version	Bp Premier version <b>Spectra SP1 Revision 1</b> 1.13.1.1078.
Notes last updated	27 May 2025
Which version can I upgrade from?	You can upgrade to version <b>Spectra SP1 Revision 1</b> (1.13.1.1078) from Bp Premier version <b>Orchid SP1 Revision 11.12.1.1023</b> or later. To check the current program version that you have installed, select <b>Help &gt; About</b> .
Which data update do I need?	You must install the <b>April 2025</b> Data Update or later before you can upgrade to version 1.13.1.1078. However, Best Practice Software recommend making sure you have the most recent Data Update installed before you run a program upgrade.
Which database and operating system versions are supported?	Windows operating system Windows 10 Pro or higher is supported. Microsoft Windows Server 2016 and higher is supported. Microsoft SQL Server 2016 and higher is supported.
When should I upgrade?	Best Practice Software recommend running the program upgrade outside of business hours.
How do I upgrade?	Brief upgrade instructions are included in this document.

---

### Copyright Statement

This material is classified as commercial-in-confidence. Unauthorised distribution of this information may constitute a breach of our Code of Conduct, and may infringe our intellectual property rights. This information is collected and managed in accordance with our [Privacy Policy](#), available on our website. © Copyright 2025.

## Upgrade to Spectra SP1 Revision 1 1.13.1.1078


---

You must have the **April 2025** data update or later installed before you upgrade to **Spectra SP1 Revision 1 1.13.1.1078**.

1. Back up your Bp Premier database before running the upgrade.
2. Synchronise all workstations in remote database mode with the Bp Premier server before you upgrade.
3. [Download and run](#) the latest Drug Update on your Bp Premier server computer.
4. Log all users out of the server and all clients before you upgrade.
5. Download and run the **Spectra SP1 Revision 1 1.13.1.1078** update on your server and all clients.
6. Log back in to Bp Premier and resume use.

## Where do I find more information?

---

Select  from the main screen of Bp Premier or select **Help > Online** from Bp Premier to open the Knowledge Base from any main screen toolbar.

## New Features

---

### Date of Birth requirements on all Prescriptions in the ACT

The [Australian Capital Territory \(ACT\)](#) has amended the [Medicines, Poisons, and Therapeutic Goods Act 2008](#) to mandate that the patient's **Date of Birth** appear on all **Prescriptions** issued, dispensed or administered in ACT from **13 June 2025**.

In Spectra SP1 Revision 1, when a prescriber logs into a **Location** with a **Postcode** in the **Australian Capital Territory (ACT)**, the patient's **Date of Birth** will be displayed on all **Prescriptions**, including those that are printed.

## Other enhancements

Bp Function	Release Notes	Key
Current Rx	<p>The <a href="#">Australian Capital Territory (ACT)</a> has amended the <a href="#">Medicines, Poisons, and Therapeutic Goods Act 2008</a> to mandate that the patient's <b>Date of Birth</b> appear on all <b>Prescriptions</b> issued, dispensed or administered in ACT from <b>13 June 2025</b>.</p> <p>When a prescriber logs into a <b>Location</b> with a <b>Postcode</b> in the <b>Australian Capital Territory (ACT)</b>, the patient's <b>Date of Birth</b> will be displayed on all <b>Prescriptions</b>, including those that are printed.</p>	79864
Bp Comms	The <b>ReceiptID</b> value for <b>Bp Comms</b> SMS delivery reports and replies has been increased.	79869

## Bug fixes

Bp Function	Release Notes	Key
Suspect Databases	Resolved an issue where if a <b>Suspect Database</b> error occurred, and the <b>Repair</b> option was selected, a connection to the server error notification was displayed, and the Repair could not be completed.	79665
Bank Accounts	Resolved an issue that prevented duplicate Provider/Practice <b>Bank Accounts</b> from appearing in the <b>Setup &gt; Configuration &gt; Bank Accounts</b> window when they had the same <b>Account Number</b> .	79994
Database Connections	Resolved an issue where if a call stack error occurs and is not closed, and a network connectivity error occurs, the application will not regain full connectivity with the server.	80075
Database Connections	Resolved an issue in which users were not warned when a network connectivity error occurred while accessing or using the <b>Bp Premier Utilities</b> .	80078

## Known issues

---

Bp Function	Issue	Fixed In	Key
Database Connections	<p>When a connection issue occurs, the <b>Retry prompt</b> appears prominently in front of the active window, preventing users from switching to other windows or using Bp Premier in the background. However, if the user opts to navigate to other windows using <b>Shortcut keys</b>, there is a risk of data loss.</p> <p><b>Workaround:</b></p> <p>To minimise the risk of data loss, users should avoid using <b>Shortcut keys</b> while Bp Premier is attempting to reconnect to the Bp Premier database server.</p>	Not yet resolved.	80048